



Active Vibration Isolation Platform SW800-4.0

- dimensions: 800x800x130mm
- weight: 105kg
- maximum load capacity: 600kg
- material of housing: aluminum
- surface: anodized
- configuration: 4 internal isolation elements (16 sensor/actuator pairs)
- electronics: 4-channel controller
- cables: 4 D-sub cables (3m long)

D-sub connectors require a space of 100mm

position of D-sub connectors

4-channel controller

		Allgemein-toleranzen ISO 2768-f		Werkstück-kanten DIN ISO 13715		Werkstoff:		Maßstab:	
		Datum		Name		Active Vibration Isolation Platform SW800-4.0 and Controller			
Gezeichnet		06.04.2020		U. Stöbener					
Kontrolliert									
Norm						Artikel-Nr.:		Zeichnungs-Nr.:	
Status		Änderungen		Datum		Name		F:\CAD\Technische Zeichnungen\MOD-Sandwich Zeichnungen\Sondermaß 800x800\SW800_4_Controller A3.idw	

A3